

Market Leading Clinical Chemistry Reagents from Europe

Comprehensive Test Menu, Calibrators & Controls

22 Test Panels, 67 Parameters

Longest On-Board and Calibration Stabilities


CORMAY

 TechnoMed

Clinical Chemistry Technology from Europe

Routine Screening

ALAT, Albumin, ALP, Amylase, ASAT, Bilirubin Total, Bilirubin Direct, Calcium, Chloride ISE, Cholesterol, HDL Direct, LDL Direct, CK, Creatinine, GGT, Glucose, Iron, LDH, Lipase, Magnesium, Phosphorus, Potassium ISE, Sodium ISE, Triglycerides, Total Protein, Urea, Uric Acid, Urine Proteins

Allergic

Total IgE

Anemia

Ferritin, Iron, Transferrin, UIBC

Bleeding disorders

Antithrombin III

Bone metabolism

ALP, Calcium, Phosphorus

Cardiac

ALAT, Apo B, ASAT, Cholesterol, HDL Direct, LDL, Direct, CK, CRP, Glucose, HBDH, Lactate, LDH, Lipoprotein (a), Myoglobin, Potassium ISE, Sodium ISE, Triglycerides

Diabetic

Cholesterol, HDL Direct, LDL Direct, Creatinine, Cystatin C, Glucose, HbA1c Direct, Lactate, Microalbumin, Triglycerides

Electrolytes

Calcium, Chloride ISE, Iron, Magnesium, Phosphorus, Potassium ISE, Sodium ISE

Hepatic

AAT, ALAT, Albumin, ALP, ASAT, Antithrombin III, Bilirubin Direct, Bilirubin Total, Ceruloplasmin, Cholinesterase, GGT, Lactate, Total Bile Acids, Total Protein, Transferrin

Infectious

AAT, AGP, ASO, Ceruloplasmin, C3, C4, CRP, Haptoglobin

Inflammatory

AAT, AGP, ASO, Ceruloplasmin, C3, C4, CRP, Ferritin, Haptoglobin, Procalcitonin, RF

Immune

C3, C4, IgA, IgG, IgM, Total IgE

Lipids

Apo B, Cholesterol, HDL Direct, LDL Direct, Lipoprotein (a), Triglycerides

Nutritional

Albumin, Iron, Magnesium, Total Protein, Transferrin, Uric Acid, Urea

Oncological

AGP, Calcium, Cholinesterase, Ferritin, IgA, IgG, IgM, Lactate, LDH, Lipase, Total Bile Acids, Total IgE, Uric Acid

Pancreatic

Amylase, Glucose, Lipase

Pediatric Screening

Albumin, AAT, ASO, Bilirubin Direct, Bilirubin Total, CRP, Creatinine, Cystatin C, IgA, IgG, IgM, Iron, Lactate, Magnesium, Phosphorus, Total IgE, Total Protein, Transferrin

Pregnancy

ALAT, ALP, ASAT, Bilirubin Total, Calcium, Creatinine, CRP, Ferritin, Glucose, GGT, Iron, Magnesium, Potassium ISE, Sodium ISE, Total Bile Acids, Transferrin, UIBC

Renal

Albumin, Calcium, Ceruloplasmin, Creatinine, Cystatin C, Glucose, Lactate, Microalbumin, Phosphorus, Potassium ISE, Sodium ISE, Urea, Uric Acid, Urine Proteins

Specific Proteins

AAT, AGP, Apo B, ASO, Antithrombin III, Ceruloplasmin, C3, C4, CRP, Cystatin C, Ferritin, Haptoglobin, HbA1c Direct, IgA, IgG, IgM, Lipoprotein (a), Microalbumin, Myoglobin, Procalcitonin, RF, Total IgE, Transferrin

Toxicological

Ethanol

COVID-19

ALAT, Albumin, ASAT, Bilirubin Total, Creatinine, CRP, Ferritin, LDH, Procalcitonin, Urea



Accent Patented Intelligence
Clinical Chemistry Diagnostic Test Results You Can Trust

NEW ACCENT SERIES ANALYZER

together with **the highest quality dedicated reagent line**, calibrators and controls provide **COMPLETE SOLUTION** suitable for small to mid-size laboratories. Fully automated, bench-top, compact size, small reagent consumption and reliable technology, are basic characteristics of the new generation of ACCENT family.

○ Accent S120

- Up to 120 T/H & 400 T/H with ISE
- 40 Disposable Cuvettes
- 8 Wavelengths

○ Accent MC240

- 240 T/H Constant & 400 T/H with ISE
- 80 Reusable Cuvettes
- Clog Detection
- Grating System with 12 Wavelengths

○ Accent M320

- Up to 320 T/H & 540 T/H with ISE
- 80 Reusable Cuvettes
- Clog Detection
- Grating System with 12 Wavelengths



Uninterrupted Run

Intelligent Probe:

- collision detection
- clog detection*
- liquid level detection
- inventory monitoring
- reagent pre-warming

Accent Intelligent Multi-Patented Leading-Edge Technology



Cost Savings

- Dedicated Applications Optimized for Maximum Performance
- Reduced Reagents Consumption – Only 100 µl Minimum Reaction Volume
- Longest Reagents On-Board and Calibration Stability (up to 12 Weeks)



Uncompromised Patient Diagnosis

- Closed & Validated Complete Clinical Chemistry Diagnostic System (Analyzer + Reagents + Calibrators + Controls)
- Patented Optical System* – Advanced Concave Grating Light-Splitting Technology*
- Waterfall* Probe Cleaning $\leq 0.05\%$ Carry-Over
- Optimized Washing Station ≤ 1 ppm Carry-Over *
- Independent Mixing Bar



Advanced Software

- Intuitive User Friendly Interface
- Intelligent Reflex Function for Conclusive Diagnosis
- Serum Index for Analyzing the Interferents (Hemolysis, Icterus and Lipemia)
- Intelligent Self Diagnostics
- Flexible During Run Reagent Loading

*available for Accent MC240 and M320