

# URIT



# TechnoMed

Compact design and easy operation

## URIT-50 Urine Analyzer



Global Diagnostics Supplier

# URIT-50

## Urine Analyzer

FDA 510 ( K ) No.K082811 Approved



### Specifications

<b>Parameters:</b> pH, Nitrite, Specific Gravity, Ascorbic Acid, Blood, Glucose, Bilirubin, Urobilinogen, Ketone, Leucocytes, Protein, Creatinine, Calcium, Microalbumin
Measurement Principle: Dual wavelength reflectance method
<b>Measuring Wavelength:</b> 3 LED Wavelength (550,620,720nm)
<b>Sample supply:</b> manual dipping
<b>Reaction time:</b> 60 seconds
<b>Throughput:</b> Normal mode, 60 samples /hour Continuous mode, 125 samples /hour
<b>Display:</b> LCD screen
<b>Memory capacity:</b> 5000 results
<b>Strip:</b> URIT 10G/11G/14G urine reagent strip.
<b>Interface:</b> Serial port(RS232), parallel port (For external printer), PS/2(For barcode reader)
<b>Work condition:</b> Environmental Temperature: 15–30°C, RH: 20%–80%
<b>Power supply:</b> Adopter 12V DC, 4.0A
<b>Dimension:</b> 290mm x 200mm x 100mm
<b>Weight:</b> 2kg



URIT URINE CONTROL

### Features

- ✓ Cold light source facility with long life and high stability
- ✓ Internal thermal printer
- ✓ Easy operation, cost-efficiency and high quality
- ✓ Rs232 serial interface to host PC/LIS, Barcode reader optional
- ✓ URIT UQ-10, UQ-11 & UQ-14 Urinalysis Control



### URIT 10G、11G、14G urine reagent strip

14G with Micro albumin, Creatinine and Calcium.

## URIT Medical Electronic Group Co., Ltd.

Head Office

Add.: No. D-07 Information Industry District,  
High-tech Zone, Guilin, Guangxi 541004,  
P. R. China

Tel: +86-773-2288586 Fax: +86-773-2288560

Shenzhen Office

Add.: Room E.J, 24th Floor, New Baohui Building, Nanhai Road,  
Nanshan District, Shenzhen, P.R. China

Tel: +86-755-26050853 Fax: +86-755-86036150  
E-mail: export@uritest.com service@uritest.com Http://www.urit.com



ISO9001  
ISO13485