

URIT

 TechnoMed

High reliability & automatization with convenience

URIT-500B Urine Analyzer



Global Diagnostics Supplier

URIT-500B Urine Analyzer

FDA 510 (K) No.K082811 Approved



Specifications

Parameters: Leukocyte, Ketone, Nitrite, Urobilinogen, Bilirubin, Protein, Glucose, Specific Gravity, Blood, pH, Ascorbic Acid, Microalbumin, Creatinine, Calcium

Measurement Principle: Dual wavelength reflectance photometry

Reaction time: 65 seconds

Throughput: 520 samples /hour

Display: LCD touch screen

Memory Capacity: 9999 results

Strip: URIT 10G/11G/14G urine reagent strip.

Interface: Serial port(RS232), parallel port (For external printer), PS/2(For barcode reader),USB Port(For data out put)

Work condition: Environmental Temperature: 15-30°C,
RH: 20%~80%,Lumination:=2000 lx

Power supply: 100~240V AC,50/60Hz

Dimension: 390mm x 340mm x 290mm(Depth×Width×Height)

Weight: 6.5kg

Features & Functions

- ▣ Friendly operation with touch screen
- ▣ Easy maintenance design with auto waste collection
- ▣ High test speed, ideal for extensive sample laboratory
- ▣ Fully automatic operation for ready-to-test
- ▣ Internal calibration for sensitivity of strips
- ▣ Rs232 serial interface to host PC/LIS, Barcode reader optional
- ▣ URIT UQ-10,UQ-11 and UQ-14 Urinalysis Control



URIT 10G、11G、14G urine reagent strip

14G with Micro albumin, Creatinine and Calcium.

URIT Medical Electronic Co., Ltd.

Head Office

Add.: No. D-07 Information Industry District,
High-tech Zone, Guilin, Guangxi 541004,
P. R. China

Tel: +86-773-2288586 Fax: +86-773-2288560

Shenzhen Office

Add.: Room E.J, 24th Floor, New Baohui Building, Nanhai Road,
Nanshan District, Shenzhen, P.R. China

Tel: +86-755-26050853 Fax: +86-755-86036150

E-mail: export@urit.com service@urit.com Http://www.urit.com

